

# The Braces KTJ®

## Joint with variable center of rotation



### Aperta

### For Functional Recovery



#### Description

During the second phase of rehabilitation, after:

- Ligament traumas
- Capsular traumas
- Meniscus traumas

#### Joint

Immobilization:

0°, 15°, 30°, 45°, 60°

Arc of controlled movement:

0°-5°-30°-45°-60°

Freedom of movement:

0° - 135°



#### Characteristics

- Permanent bacteriostatic action
- Hypo-allergenic and anti-odour
- Breathes with the skin
- Absorbs perspiration
- Dries quickly
- Dirt resistant and washable
- Doesn't stick
- Anti-static action

#### Size

XS: 33 - 38 cm.

S: 39 - 44 cm.

M: 45 - 51 cm.

L: 52 - 56 cm.

XL: more than 57 cm.

Right and Left



The circumference of the thigh should be measured 15 cm. from the top of the kneecap.