

The Braces KTJ®

Joint with variable center of rotation



Atleta

Polyvalent Brace



Description

During the second phase of rehabilitation, after:

- Ligament traumas
- Capsular traumas
- Meniscus traumas

For resumption sport activity

Joint

Immobilization:

0°, 15°, 30°, 45°, 60°

Arc of controlled movement:

0°-5°-30°-45°-60°

Freedom of movement:

0° - 135°



Characteristics

- Permanent bacteriostatic action
- Hypo-allergenic and anti-odour
- Breathes with the skin
- Absorbs perspiration
- Dries quickly
- Dirt resistant and washable
- Doesn't stick
- Anti-static action

Size

One Size: 40 - 54 cm.

Right and Left